

TSS Test

A Citizen Scientist Adaptation

Test Procedure

- Collect water sample
- Pre-dry and pre-weigh filter paper
- Filter a known volume of the sample
(if from large sample, pull aliquot while agitating)
- Dry filter paper to constant weight
- Weigh using analytical balance
- Calculate TSS as mg of solids per liter of sample

TSS Test

- The test is straightforward.
- However, if useable results are to be obtained good lab technique is essential.
- A simple lab setup can be put together at a modest cost to give results approaching those from a professional lab.
- Parts which follow can be purchased on Amazon unless noted

Measure out aliquot of sample

Agitate sample well before measuring out an aliquot- **immediately**
May want different size graduated cylinders



by [EISCO](#)

Graduated Cylinder, 100ml -
Polypropylene, Class B - Octagonal
Base, Autoclavable - Industrial Quality,
US Sourced Materials - Eisco Labs

★★★★☆ [46 customer reviews](#)

| [7 answered questions](#)

Price: **\$8.49** ✓ [prime](#) FREE One-Day

New (2) from **\$8.49** ✓ [prime](#)

Specifications for this item

Filter Sample- Option A

- Follow as closely as possible standard TSS test protocol
- Filter paper of the glass fiber/non-binder type
- Filter paper should be pre-washed with de-ionized or distilled water
- Washed filter paper must be pre-dried till no further weight loss occurs before determining starting weight
- Equipment:
 - Filter funnel
 - Filter flask
 - Filter paper
 - Vacuum source

Filter Kit



by D DOLITY

D DOLITY 1000ml Filtration Apparatus
Vacuum Lab Filtering Unit Flask Funnel
and Clamp, Compatible with 42-60mm
Membrane Filter Paper

[Be the first to review this item](#)

Price: **\$35.28** & **FREE Shipping**

New (1) from \$35.28 & FREE shipping.

Specifications for this item



by [Whatman](#)

**Whatman 1821-047 Glass Microfiber
Binder Free Filter, 1 Micron, 12
s/100mL Flow Rate, Grade GF/B, 4.7cm
Diameter (Pack of 100)**

[Be the first to review this item](#)

Price: **\$104.75** (\$1.05 / Filters) & **FREE Shipping**. [Details](#)

New (3) from \$98.00 & FREE shipping.

Specifications for this item



Simsii PTFE Membrane Filter, Hydrophobic, Diameter 47 mm, Pore Size 0.45 μm , Pack of 200

by Simsii

[Be the first to review this item](#)

Price: **\$129.99** (\$6.50 / 10 Items) & **FREE Shipping**

- Hydrophobic PTFE
 - High chemical and temperature resistance
 - Resistivity against all know solvents
-

Specifications for this item

Filter Paper

Example from Hach

TSS Glass Fiber Filter, Pore Size 1.5 μm , Diameter 47 mm, 100/pk

- Overview
- Details
- Parameter/Reagent
- Downloads
- Accessories
- Similar Products




» Gallery

Product #: 253000

USD Price: \$45.89

Available

Quantity



One Time 

Add to Cart 



Very low cost Vacuum pump About \$75 on Amazon



Kozyvacu TA350 Single-Stage Rotary Vane Vacuum Pump for HVAC/Auto AC Refrigerant Recharging, Wine Degassing, Milking, Medical, Food Processing

Size:KZTA350001



Note- should provide a vacuum break

Filter Sample- Option B

- Recommend as fine filter paper as possible to get “reasonable” filtering times. Bill Ward has used 10 micron pore size.
- Filter paper MUST be pre-washed with de-ionized or distilled water to remove all water soluble materials
- Washed filter paper must be pre-dried till no further weight loss occurs before determining starting weight
- Equipment:
 - Buchner funnel
 - Filter flask
 - Vacuum source
 - Filter paper

Filter Flask and Buchner Funnel



Cost about \$60

For safety tape vacuum flask with electrical tape in case of implosion



by [Whatman](#)

Whatman 1003-070 Quantitative Filter Paper Circles, 6 Micron, 26 s/100mL/sq inch Flow Rate, Grade 3, 70mm Diameter (Pack of 100)

[Be the first to review this item](#)

Was: ~~\$31.69~~

Price: **\$23.79** (\$0.24 / Filter Paper Circles) ✓ **prime**

FREE Shipping on orders over \$25—or get **FREE**

Two-Day Shipping with [Amazon Prime](#)

You Save: **\$7.90 (25%)**

New (2) from \$13.92 + \$6.29 shipping

Super low cost vacuum source About \$25 on Amazon



Wash bottle for DI/Distilled water
About \$5 on Amazon



Tweezers about \$5-10 on Amazon

Never touch filter paper with
fingers!!



Drying-Low cost air oven (About \$500 on Amazon)



Determine conditions to dry filter/sample to constant weight at 105 C
Aluminum weighing dish/boat recommended (also on Amazon)

Super Low Cost Drying Equipment

- Oven in a kitchen
- Infrared heat lamp

Must determine conditions to dry filter/sample to constant weight

Dessicator for cooling dried samples and storage of prewashed and dried filters



Bel-Art "Space Saver"
Polycarbonate Vacuum
Desiccator with Clear
Polycarbonate Bottom; 0.31
cu. ft. (F42027-0000)

Size: 0.31 cu. ft.

Material Type: Polycarbonate with All-Clear
Polycarbonate Bottom



\$137 at Amazon



DRY&DRY Half(1/2) Gallon Replacement Desiccant Indicating Silica Gel Beads Reusable, Blue

★★★★★ 26

\$22⁹⁹

✓prime Get it as soon as **Tomorrow,**

Jun 5

FREE Shipping by Amazon



Disposable Aluminum 70 mm Smooth Wall Weighing Dish- Pack of 100

★★★★☆ ~ 3

\$22⁰⁰ (\$2.20/10 Items)

✓prime Get it as soon as **Thu, Jun 6**

FREE Shipping on orders over \$25
shipped by Amazon

Only 19 left in stock - order soon.

Analytical Balance weigh to 0.1 mg (About \$500 on Amazon)



Sponsored ⓘ

120g 0.0001g 0.1mg Digital Analytical Balance Scale for Laboratories from U.S. Solid

★★★★★ ∨ 2

\$495¹¹

✓prime FREE delivery

TSS Calculation

Dry suspended solids weight = filter final weight – filter starting weight
(convert to milligrams)

Convert sample volume in cc or ml to liters (divide by 1000)

TSS = solids weight/ liters of sample: units of mg/liter

Stroud Filtration membrane filter 0.45 micron Millipore					10 micron GFF filter			
Stroud Filter #	ID	Vol Filtered (mL)	total mass sediment (mg)	mass corrected for balance drift	TSS (mg/L)	TSS (mg/L) (info from Bill Ward)	Volume Filtered (info from Bill Ward)	(10 um filter weight) / (0.45um filter weight)
M8315		255	19.85	20.00	78.43	56		0.714
M8324		294	13.93	13.58	46.19	33		0.71443299
M8325		122	92.34	92.98	762.13	654		0.858120026